



Sterilization is more effective

# K9 Automatic Temperature Measurement & Disinfection Machine

Safer without contact



Large-capacity



Automatic dispensing



Temperature measurement



Alarm Function



High precision



°C/°F






# Automatic temperature measurement and disinfection

## Effectively avoid cross infection

Align the temperature measuring probe and the liquid outlet with the palm of your hand to detect the temperature and automatically sterilize the liquid



Extremely fast temperature measurement



Infrared sensor

# Automatic alarm function

When the measured body temperature exceeds the preset alarm temperature, the device will issue an alarm to remind the surrounding people



High temperature alarm



# Accurate measurement

accurate measurement, fast and efficient, stronger  
environmental adaptability



# Large-capacity disinfectant used longer

1000ML large amount of disinfectant solution, effectively inhibit bacteria, the inhibition rate is as high as 99.9%.  
Effectively take care of you and your family's health



# °C/°F temperature unit can be selected

Switch easily






# Self-developed core motherboard chip

High-performance, low-power self-developed motherboard  
greatly improves product stability and operating capability

The inner core is more powerful

K9



# Smart power saving & auto wake up

K9 has intelligent power saving and automatic wake-up function, which can save more resources

Low energy consumption, more economical



# Visualization design of disinfectant

The side adopts a transparent design to visualize the use of disinfectant, and add liquid at any time according to the use





**1. 條款**

1. 所有權保留：在收到全部付款前，所有權保留在賣方。
2. 付款：請在收到貨物後70天內付款。
3. 運費：運費由買方自理。
4. 保險：保險由買方自理。

**2. 付款方式**


1. Western Union / T / PayPal
2. 30天內付款
3. 30天內付款

**3. 運輸**

運輸方式：/ / / /  
TNT / FedEx / DHL / UPS / EMS  
1-3天  
5-7天  
OEM 7-10天

**4. 保固**

1. Warranty 12個月
2. 保固條款
3. 保固條款



# Technical Data



Accuracy	+/-0.2°C (34~45°C, put it in the operating environment for 30 minutes before use)		
Alarm type	Flashing+ bleeping “DIDIDI”	Input	DC 5V 1A
Measuring distance	5cm~10cm	Product power	3W
LED display	Number	Liquid type	Hand sanitizer
Interface	USB Type-C Port	<p style="text-align: center; color: #008080; font-weight: bold;">Advantages</p> <ul style="list-style-type: none"> <li>• Automatic induction no-clean disinfection soap dispenser, more convenient to use</li> <li>• Large digital display, high-definition display</li> <li>• High-precision sensor, stable and reliable performance</li> <li>• With intelligent point saving and automatic wake-up function</li> <li>• Strong adaptability to environment temperature</li> <li>• Large-capacity soap liquid, longer life cycle</li> </ul>	
Bottle Capacity	1000mL		
Mounting Type	Wall-hanging, Tripod-fixed		
Working Temperature	8.10°C~40°C ( 15°C~35°C Recommended)		
Measuring range	0~50°C		
Standby time	5s for non-working, dynamic loop with white dot		
Dimension	119x133x280mm		
Net weight	702g		

# Box Contents



K9



User Guide




USB cable&charger\*1



Screw&rubber plug\*2

# Appearance





# Strong environmental adaptability



Hospital



School



Community



Factory



Supermarket



Station



Hotel



Exhibition





Market



□□□□	1□100	101□5000	5001□10000	> 10000
EST□□□□□□□□□□	5	7	10	□□□□

□□□





Q: 请问这款皂液器的容量是多少？  
A: 这款皂液器的容量是3000毫升。

Q: 请问这款皂液器的使用寿命有多长？  
A: 这款皂液器的使用寿命可达24个月。

Q

A EXW FOB CIF

Q

A 2 OEM / OEM