

00 0000 000 00 00 00 00,Bluetooth 00 00 00 00 LED Metene 500dl FDA 000 0 00 00 00 00 000000

00:

6 00000 00 00, 00 0 00000 00, 00 00 OLED 000000, 00 0 00 00

00 0000, 00 00 : 0.2 %

000 00000 00 000 0000000.

00 000 00 000, 000000 00 00000 00 00, 00 2 00 AAA 1.5V 00000 000, 20 00 00 00 00 00, 00 000 000000.

000:

00 00000 000 00 00 000-SpO2	
00 00:	35 ~ 100 %
00:	1%
00 000 :	2 % (70-100 %),00000 00 (≤70 %)
00 00 :	50 % -100 %
0000	
00 00:	30-250BPM
00:	1BPM
00 000 :	2BPM
000:	FDA



SpO<sub>2</sub>%

97

PR bpm

65

PI%

14.5



**FINGERTIP PULSE OXIMETER**

- SIMPLE OPERATION
- DISPLAY FAST
- EASY TO CARRY
- HIGH ACCURACY

□□ □□ □□ □□ □□ □□ □□

□□ □ □□

□□ (L × W × H)	58mm × 37mm × 32mm
□□□	54g
□□ □□ □□ (L × W × H)	86mm × 60mm × 38mm
□□ □□ □□ (L × W × H)	430mm × 340mm × 200mm
□□□ (100pcs)	7Kg (□□□ □□)
□□	□□

## Shipping service



□□ : □□□ □□ □□ **TT** □□□ **30 %**, □□□ □□ **TT** □□□ **70 %** □□

□□ : **30 %** □□□□□ □□ **20days** □□. **OEM**

□□ □□□□ **60** □□□□.