

Productomschrijving

UV-licht Handheld UV-sterilisator met 18 stuks importeren led voor kamer desinfectie

Productnaam	UV-sterilisator	Aantal lampkralen	18 LED
Materiaal	buikspieren	Toepassing	Desinfectie
Kracht	7W		Mijten verwijderen
Grootte	20,7 * 4,5 * 2,8 cm		Steriliseren
Voorzien zijn van	LED UV-desinfectie, sterilisatiegraad van 99,9%		

Verpakking en levering

Verkoopeenheden	Enkel voorwerp
Enkele pakketgrootte	20X5X3 cm
Enkel brutogewicht	0,002 kg



18pcs of Importing **LED** beads

99.9% sterilization rate for
home & travelling disinfection



Gravity sensor setting

Automatically turns off UV light

When separating the upper part

99.9%

LED UV disinfection

Destroyed over 99.9% of germs

Ultraviolet disinfection technology is commonly used in hospitals, bacteriological research, military and aerospace fields.

It uses a purely physical sterilization method and irradiates with a wavelength of 260-280nm.

Can instantly destroy DNA and RNA of microorganisms such as bacteria.

Make it lose its ability to reproduce and survive.



Light in weight, 200g

1s disassembly and assembly, portable

Chargeable battery, USB charging

Full charging for using 60mins

3000mAh lithium battery



Adjustable UV sanitizing wand


Long & Short size for using different place

Updated from tradition UV sterilizer,

let it more suitable and useful to lives

Short sterilizer is used to :
mask, cellphone, cup, corners,
all-around, more energy-saving.

Combined long sterilizer is used to:
car, bedroom, etc..big space,
experience the small size bringing
you a safety and health living



Fight Against Bacteria

Respiratory
disease

Bacteria
Breeding
ground

Bacteria
Often allergic

Sneezing



Details Afbeeldingen

Gebruik altijd en overal in je leven om jezelf een veilige en gezonde leefruimte te geven

Waar te gebruiken I:
Gebruik voor uw maskers

Slightly sweep 5 times, quick disinfection
Reuse to save resources



Waar II te gebruiken:
Gebruik voor uw mobiele telefoon



Waar III te gebruiken:
Gebruik voor uw servies

Short sterilizer is used to disinfect small stuff

Safe and save energy



Waar te gebruiken IV
Gebruik voor uw slaapkamer



Waar V te gebruiken:
Gebruik voor kinderwagen



Waar VI te gebruiken
Gebruik voor uw toilet




Productpresentatie

18 PCS UV lamp beads for
10 years life

99.9%

Strong sterilization
Best effectivity
Most safety



Certificaat



广微测
Gmicro Testing



201819C00883



中国认可
国际互认
检测
TESTING
CNAS L1747

广东省微生物分析检测中心

GUANGDONG DETECTION CENTER OF MICROBIOLOGY

分析检测结果

ANALYSIS AND TEST RESULT

报告编号 (Report No.): 2020FM06103R01D

1.样品前处理:以玻璃为载体,作用距离3cm,照射时间为5分钟。

Sample pretreatment: Took the glass as the carrier, and the action distance was 3cm, Irradiation and the action time was 5 minutes.

2.试验结果 Test Result:

试验菌株 Test germ	试验 组别 Test group	对照组平均菌落数 Average cfu of positive controls (cfu/片)	试验组平均菌落数 Average cfu of testing groups (cfu/片)	杀灭率 Killing rate (%)	杀灭对数值 Sterilization logarithm
大肠杆菌 (<i>Escherichia coli</i>) 8099	1	5.0×10^6	1.8×10^3	99.96	3.44
	2	4.7×10^6	1.8×10^3	99.96	3.41
	3	4.8×10^6	1.9×10^3	99.96	3.40
金黄色葡萄球菌 (<i>Staphylococcus aureus</i>) ATCC 6538	1	1.6×10^6	7.0×10^2	99.96	3.35
	2	1.4×10^6	7.7×10^2	99.95	3.26
	3	1.3×10^6	7.5×10^2	99.94	3.23

(以下空白 Blank below)



Welkom bij O-leading

Omdat het coronavirus zich wereldwijd verspreidde, vooral tijdens Jan tot februari in China, hadden we zelfs een groot tekort aan maskers en soorten goederen tegen COVID-19! Nu wordt China beter en het is de O-leidende verantwoordelijkheid om de wereldvrienden met deze goederen te helpen, raadpleeg alstublieft meer over onze website [Antivirus](#)



White Believe!
FIGHTING COME ON